

产品名称: Protein Jumonji Rabbit Monoclonal Antibody
产品货号: AMRe85721



产品概述 (Summary)

产品名称 (Production Name)	Protein Jumonji Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Protein Jumonji
别名 (Alternative Names)	JARD2; JARID2; JMJ
基因 ID (Gene ID)	3720.0
蛋白 ID (SwissProt ID)	Q92833.A synthetic peptide of human Jarid2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 139 kDa; Observed MW: 139 kDa

研究背景 (Background)

Regulator of histone methyltransferase complexes that plays an essential role in embryonic development, including heart


产品名称: Protein Jumonji Rabbit Monoclonal Antibody
产品货号: AMRe85721



and liver development, neural tube fusion process and hematopoiesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Jarid2 in K562, C6 lysates using Protein Jumonji antibody.

注意事项 (Note)

For research use only .