

产品名称: Kaiso Rabbit Monoclonal Antibody
产品货号: AMRe85725



产品概述 (Summary)

产品名称 (Production Name)	Kaiso Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Kaiso
别名 (Alternative Names)	Kaiso; Zbtb33; ZNF Kaiso; ZNF348
基因 ID (Gene ID)	10009.0
蛋白 ID (SwissProt ID)	Q86T24.A synthetic peptide of human Kaiso

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 74 kDa; Observed MW: 90 kDa

研究背景 (Background)

Transcriptional regulator with bimodal DNA-binding specificity. Binds to methylated CpG dinucleotides in the consensus

产品名称: Kaiso Rabbit Monoclonal Antibody
产品货号: AMRe85725




sequence 5'-CGCG-3' and also binds to the non-methylated consensus sequence 5'-CTGCNA-3'.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis of Kaiso in MCF-7, C6 lysates using Kaiso antibody.

注意事项 (Note)

For research use only .