

产品名称: KAT1 Rabbit Monoclonal Antibody
产品货号: AMRe85728



产品概述 (Summary)

产品名称 (Production Name)	KAT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KAT1
别名 (Alternative Names)	hat1; KAT1
基因 ID (Gene ID)	8520.0
蛋白 ID (SwissProt ID)	O14929.Recombinant protein of human KAT1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 45 kDa

研究背景 (Background)

Acetylates soluble but not nucleosomal histone H4 at 'Lys-5' (H4K5ac) and 'Lys-12' (H4K12ac) and, to a lesser extent,


产品名称: KAT1 Rabbit Monoclonal Antibody
产品货号: AMRe85728




acetylates histone H2A at 'Lys-5' (H2AK5ac). Has intrinsic substrate specificity that modifies lysine in recognition sequence GXGKXG.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of KAT1 in K562, C6, 3T3, HeLa lysates using KAT1 antibody.



Immunocytochemistry analysis of KAT1(green) in HeLa using KAT1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KAT1/HAT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: KAT1 Rabbit Monoclonal Antibody
产品货号: AMRe85728



注意事项 (Note)

For research use only .