

产品名称: KCTD9 Rabbit Monoclonal Antibody
产品货号: AMRe85731



产品概述 (Summary)

产品名称 (Production Name)	KCTD9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KCTD9
别名 (Alternative Names)	BTB/POZ domain-containing protein KCTD9
基因 ID (Gene ID)	54793.0
蛋白 ID (SwissProt ID)	Q7L273.A synthetic peptide of human KCTD9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

研究背景 (Background)

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex, which mediates the


产品名称: KCTD9 Rabbit Monoclonal Antibody
产品货号: AMRe85731




ubiquitination of target proteins, leading to their degradation by the proteasome.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of KCTD9 in K562, HeLa lysates using KCTD9 antibody.



Immunocytochemistry analysis of KCTD9(green) in HeLa using KCTD9 antibody, and DAPI(blue)



Immunocytochemistry analysis of KDM1/LSD1(green) in HeLa using KDM1/LSD1 antibody, and DAPI(blue)

产品名称: KCTD9 Rabbit Monoclonal Antibody
产品货号: AMRe85731



注意事项 (Note)

For research use only .