

产品名称: KGF Rabbit Monoclonal Antibody
产品货号: AMRe85734



产品概述 (Summary)

产品名称 (Production Name)	KGF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KGF
别名 (Alternative Names)	Fgf7; Fibroblast growth factor 7; HBGF7; Heparin binding growth factor 7; Keratinocyte growth factor
基因 ID (Gene ID)	2252.0
蛋白 ID (SwissProt ID)	P21781.A synthetic peptide of human KGF

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 28 kDa

研究背景 (Background)

产品名称: KGF Rabbit Monoclonal Antibody
产品货号: AMRe85734




Growth factor active on keratinocytes. Possible major paracrine effector of normal epithelial cell proliferation.

研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of KGF in rat Brain, C6, 3T3, HeLa lysates using KGF antibody.

注意事项 (Note)

For research use only .