

产品名称: LIS1 Rabbit Monoclonal Antibody
产品货号: AMRe85753



产品概述 (Summary)

产品名称 (Production Name)	LIS1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Rat, Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LIS1
别名 (Alternative Names)	MDS; LIS1; LIS2; MDCR; NudF; PAFAH
基因 ID (Gene ID)	5048.0
蛋白 ID (SwissProt ID)	P43034.A synthetic peptide of human LIS1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Required for proper activation of Rho GTPases and actin polymerization at the leading edge of locomoting cerebellar


产品名称: LIS1 Rabbit Monoclonal Antibody
产品货号: AMRe85753



neurons and postmigratory hippocampal neurons in response to calcium influx triggered via NMDA receptors.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of LIS1 in C6, CHO-K1 lysates using LIS1 antibody.

注意事项 (Note)

For research use only .