

产品名称: FABP1 Rabbit Monoclonal Antibody
产品货号: AMRe85754



产品概述 (Summary)

产品名称 (Production Name)	FABP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABP1 FABP 1; FABP1; FABP1; FABPL; FABPL_HUMAN; Fatty Acid Binding Protein 1; Fatty acid binding protein 1 liver; Fatty Acid Binding Protein; Fatty acid-binding protein 1; Fatty acid-binding protein; Fatty acid-binding protein liver; L FABP; L-FABP; liver; Liver-type fatty acid-binding protein.
别名 (Alternative Names)	
基因 ID (Gene ID)	2168.0
蛋白 ID (SwissProt ID)	P07148.Recombinant protein of human liver FABP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa; Observed MW: 14 kDa

产品名称: FABP1 Rabbit Monoclonal Antibody
产品货号: AMRe85754




研究背景 (Background)

Binds free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Liver FABP in mouse liver lysates using FABP1 antibody.

注意事项 (Note)

For research use only .