

产品名称: Lumican Rabbit Monoclonal Antibody
产品货号: AMRe85761



产品概述 (Summary)

产品名称 (Production Name)	Lumican Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lumican
别名 (Alternative Names)	LDC; LUM; Lumican; SLRR2D
基因 ID (Gene ID)	4060.0
蛋白 ID (SwissProt ID)	P51884.A synthetic peptide of human Lumican

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 60 kDa

研究背景 (Background)

Belongs to the small leucine-rich proteoglycan (SLRP) family. SLRP class II subfamily. Contains 11 LRR (leucine-rich) repeats.


产品名称: Lumican Rabbit Monoclonal Antibody
产品货号: AMRe85761




Contains 1 LRRNT domain.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Lumican in K562, C6, 3T3, HeLa lysates using Lumican antibody.



Immunocytochemistry analysis of Lumican(green) in HeLa using Lumican antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Lumican antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Lumican Rabbit Monoclonal Antibody
产品货号: AMRe85761



注意事项 (Note)

For research use only .