

产品名称: MAS1L Rabbit Monoclonal Antibody
产品货号: AMRe85769



产品概述 (Summary)

产品名称 (Production Name)	MAS1L Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAS1L
别名 (Alternative Names)	MRG; MAS-L; dJ994E9.2
基因 ID (Gene ID)	116511.0
蛋白 ID (SwissProt ID)	P35410.A synthetic peptide of human MAS1L

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 42 kDa

研究背景 (Background)

Predicted to enable G protein-coupled receptor activity. Predicted to be involved in G protein-coupled receptor signaling


产品名称: MAS1L Rabbit Monoclonal Antibody
产品货号: AMRe85769



pathway. Located in cytosol; nucleoplasm; and plasma membrane.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MAS1L in K562, rat Brain, C6, 3T3, HeLa lysates using MAS1L antibody.

注意事项 (Note)

For research use only .