

产品名称: M-CSF Receptor Rabbit Monoclonal Antibody
产品货号: AMRe85774



产品概述 (Summary)

产品名称 (Production Name)	M-CSF Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	M-CSF Receptor
别名 (Alternative Names)	M-CSF Receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD115
基因 ID (Gene ID)	12978.0
蛋白 ID (SwissProt ID)	P09581.Recombinant protein of mouse MCSF Receptor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 109 kDa; Observed MW: 109 kDa

研究背景 (Background)

Tyrosine-protein kinase that acts as cell-surface receptor for CSF1 and IL34 and plays an essential role in the regulation of


产品名称: M-CSF Receptor Rabbit Monoclonal Antibody
产品货号: AMRe85774



survival, proliferation and differentiation of hematopoietic precursor cells, especially mononuclear phagocytes, such as macrophages and monocytes. Promotes the release of proinflammatory chemokines in response to IL34 and CSF1, and thereby plays an important role in innate immunity and in inflammatory processes. Plays an important role in the regulation of osteoclast proliferation and differentiation, the regulation of bone resorption, and is required for normal bone and tooth development.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of M-CSF Receptor in K562 lysates using M-CSF Receptor antibody.

注意事项 (Note)

For research use only .