

**产品名称: Methylmalonyl Coenzyme A mutase Rabbit
Monoclonal Antibody**
产品货号: AMRe85782



产品概述 (Summary)

产品名称 (Production Name)	Methylmalonyl Coenzyme A mutase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Methylmalonyl Coenzyme A mutase
别名 (Alternative Names)	MCM
基因 ID (Gene ID)	4594.0
蛋白 ID (SwissProt ID)	P22033.A synthetic peptide of human Methylmalonyl Coenzyme A mutase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 83 kDa; Observed MW: 83 kDa

研究背景 (Background)


**产品名称: Methylmalonyl Coenzyme A mutase Rabbit
Monoclonal Antibody**
产品货号: AMRe85782




Involved in the degradation of several amino acids, odd-chain fatty acids and cholesterol via propionyl-CoA to the tricarboxylic acid cycle. MCM has different functions in other species.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Methylmalonyl Coenzyme A mutase in HeLa, rat Brain lysates using Methylmalonyl Coenzyme A mutase antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Methylmalonyl Coenzyme A mutase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .