

产品名称: MIF Rabbit Monoclonal Antibody
产品货号: AMRe85785



产品概述 (Summary)

产品名称 (Production Name)	MIF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MIF MIF; GLIF; MMIF; Macrophage migration inhibitory factor; MIF; Glycosylation-inhibiting factor; GIF; L-dopachrome isomerase; L-dopachrome tautomerase; Phenylpyruvate tautomerase
别名 (Alternative Names)	
基因 ID (Gene ID)	4282.0
蛋白 ID (SwissProt ID)	P14174.Recombinant protein of human MIF

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa; Observed MW: 12 kDa

产品名称: MIF Rabbit Monoclonal Antibody
产品货号: AMRe85785




研究背景 (Background)


MIF promotes the systemic inflammatory response by counter-regulating glucocorticoid-mediated inhibition of immune-cell activation and proinflammatory cytokine production. MIF may mediate tissue destruction through the induction of proteinases.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MIF in Jurkat lysates using MIF antibody.



Western blot analysis of MIF in C6, 3T3 lysates using MIF antibody.

注意事项 (Note)

For research use only .