

**产品名称: Mitochondrial Ribosomal Protein L11 Rabbit
Monoclonal Antibody**
产品货号: AMRe85786



产品概述 (Summary)

产品名称 (Production Name)	Mitochondrial Ribosomal Protein L11 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mitochondrial Ribosomal Protein L11
别名 (Alternative Names)	MRPL11; CGI-113; 39S ribosomal protein L11; mitochondrial; L11mt; MRP-L11
基因 ID (Gene ID)	65003.0
蛋白 ID (SwissProt ID)	Q9Y3B7.A synthetic peptide of human Mitochondrial ribosomal protein L11

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)


**产品名称: Mitochondrial Ribosomal Protein L11 Rabbit
Monoclonal Antibody**
产品货号: AMRe85786



This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Mitochondrial ribosomal protein L11 in HeLa lysates using Mitochondrial Ribosomal Protein L11 antibody.

注意事项 (Note)

For research use only .