

产品名称: MKLP1 Rabbit Monoclonal Antibody
产品货号: AMRe85787



产品概述 (Summary)

产品名称 (Production Name)	MKLP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MKLP1
别名 (Alternative Names)	CHO1; KNSL5; MKLP1; MKLP-1
基因 ID (Gene ID)	9493.0
蛋白 ID (SwissProt ID)	Q02241.A synthetic peptide of human MKLP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 110 kDa; Observed MW: 110 kDa

研究背景 (Background)

Plus-end-directed motor enzyme that moves antiparallel microtubules in vitro. Localizes to the interzone of mitotic


产品名称: MKLP1 Rabbit Monoclonal Antibody
产品货号: AMRe85787



spindles.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MKLP1 in HeLa lysates using MKLP1 antibody.

注意事项 (Note)

For research use only .