

产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe85790



产品概述 (Summary)

产品名称 (Production Name)	MLKL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLKL
别名 (Alternative Names)	MLKL
基因 ID (Gene ID)	197259.0
蛋白 ID (SwissProt ID)	Q8NB16.Recombinant protein of human MLKL

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

研究背景 (Background)

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following


产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe85790




phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage. Does not have protein kinase activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MLKL in HT-29 lysates using MLKL antibody.



Western blot analysis of MLKL in U251MG, HeLa lysates using MLKL antibody.

注意事项 (Note)

For research use only .