

产品名称: MMP9 Rabbit Monoclonal Antibody
产品货号: AMRe85795



产品概述 (Summary)

产品名称 (Production Name)	MMP9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MMP9
别名 (Alternative Names)	92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB
基因 ID (Gene ID)	17395.0
蛋白 ID (SwissProt ID)	P41245.Recombinant protein of mouse MMP9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 92 kDa

研究背景 (Background)

Matrix metalloproteinase that plays an essential role in local proteolysis of the extracellular matrix and in leukocyte


产品名称: MMP9 Rabbit Monoclonal Antibody
产品货号: AMRe85795



migration (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MMP9 in rat lung lysates using MMP9 antibody.

注意事项 (Note)

For research use only .