

产品名称: MRPS26 Rabbit Monoclonal Antibody
产品货号: AMRe85799



产品概述 (Summary)

产品名称 (Production Name)	MRPS26 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPS26
别名 (Alternative Names)	GI008; MRPS13; RPMS13; MRP-S13; MRP-S26; NY-BR-87; C20orf193; dJ534B8.3
基因 ID (Gene ID)	64949.0
蛋白 ID (SwissProt ID)	Q9BYN8.Recombinant protein of human MRPS26

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

研究背景 (Background)


产品名称: MRPS26 Rabbit Monoclonal Antibody
产品货号: AMRe85799




Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MRPS26 in K562, HeLa lysates using MRPS26 antibody.



Immunocytochemistry analysis of MRPS26(green) in HeLa using MRPS26 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .