

产品名称: MS4A14 Rabbit Monoclonal Antibody
产品货号: AMRe85802



产品概述 (Summary)

产品名称 (Production Name)	MS4A14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MS4A14
别名 (Alternative Names)	Testis development protein NYD-SP21
基因 ID (Gene ID)	84689.0
蛋白 ID (SwissProt ID)	Q96JA4.Recombinant protein of human MS4A14

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa

研究背景 (Background)


Membrane-spanning 4-domains subfamily A member 14

产品名称: MS4A14 Rabbit Monoclonal Antibody
产品货号: AMRe85802



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MS4A14 in Jurkat, rat Brain, C6 lysates using MS4A14 antibody.

注意事项 (Note)

For research use only .