

产品名称: MST4 Rabbit Monoclonal Antibody
产品货号: AMRe85805



产品概述 (Summary)

产品名称 (Production Name)	MST4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MST4
别名 (Alternative Names)	MASK; MST4
基因 ID (Gene ID)	51765.0
蛋白 ID (SwissProt ID)	Q9P289.A synthetic peptide of human MST4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 52 kDa

研究背景 (Background)

Mediator of cell growth (PubMed:11641781, PubMed:17360971). Modulates apoptosis (PubMed:11641781,

产品名称: MST4 Rabbit Monoclonal Antibody
产品货号: AMRe85805




PubMed:17360971). In association with STK24 negatively regulates Golgi reorientation in polarized cell migration upon RHO activation (PubMed:27807006).


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of MST4 in K562, C6 lysates using MST4 antibody.



Immunocytochemistry analysis of MST4(green) in HeLa using MST4 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .