

产品名称: Mucin 4 Rabbit Monoclonal Antibody
产品货号: AMRe85811



产品概述 (Summary)

产品名称 (Production Name)	Mucin 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mucin 4
别名 (Alternative Names)	Mucin-4; Ascites sialoglycoprotein; ASGP; Pancreatic adenocarcinoma mucin; Testis mucin; Tracheobronchial mucin
基因 ID (Gene ID)	4585.0
蛋白 ID (SwissProt ID)	Q99102.Recombinant protein of human MUC4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 542 kDa; Observed MW: 120 kDa

研究背景 (Background)


产品名称: Mucin 4 Rabbit Monoclonal Antibody
产品货号: AMRe85811



May play a role in tumor progression. Ability to promote tumor growth may be mainly due to repression of apoptosis as opposed to proliferation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MUC4 in PC-3 lysates using Mucin 4 antibody.

注意事项 (Note)

For research use only .