

产品名称: MVP Rabbit Monoclonal Antibody
产品货号: AMRe85814



产品概述 (Summary)

产品名称 (Production Name)	MVP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MVP
别名 (Alternative Names)	LRP; MVP; VAULT1
基因 ID (Gene ID)	9961.0
蛋白 ID (SwissProt ID)	Q14764.A synthetic peptide of human MVP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 110 kDa

研究背景 (Background)

Required for normal vault structure. Vaults are multi-subunit structures that may act as scaffolds for proteins involved in


产品名称: MVP Rabbit Monoclonal Antibody
产品货号: AMRe85814



signal transduction. Vaults may also play a role in nucleo-cytoplasmic transport.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MVP in U87-MG lysates using MVP antibody.

注意事项 (Note)

For research use only .