

产品名称: Myeloperoxidase Rabbit Monoclonal Antibody
产品货号: AMRe85819



产品概述 (Summary)

产品名称 (Production Name)	Myeloperoxidase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Myeloperoxidase
别名 (Alternative Names)	Myeloperoxidase; MPO
基因 ID (Gene ID)	4353.0
蛋白 ID (SwissProt ID)	P05164.A synthetic peptide of human Myeloperoxidase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 80-90 kDa

研究背景 (Background)

Myeloperoxidase (MPO) is a peroxidase enzyme that is part of the host defense system of polymorphonuclear leukocytes.


产品名称: Myeloperoxidase Rabbit Monoclonal Antibody
产品货号: AMRe85819




Heightened MPO levels have been associated with tissue damage and a number of pathological conditions.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Myeloperoxidase in 3T3 lysates using Myeloperoxidase antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Myeloperoxidase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .