

产品名称: NDUFA1 Rabbit Monoclonal Antibody
产品货号: AMRe85830



产品概述 (Summary)

产品名称 (Production Name)	NDUFA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFA1
别名 (Alternative Names)	MWFE; ZNF183; CI-MWFE
基因 ID (Gene ID)	4694.0
蛋白 ID (SwissProt ID)	O15239.A synthetic peptide of human NDUFA1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 8 kDa; Observed MW: 8 kDa

研究背景 (Background)

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not


产品名称: NDUFA1 Rabbit Monoclonal Antibody
产品货号: AMRe85830




to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NDUFA1 in K562, HeLa lysates using NDUFA1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using NDUFA1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .