

产品名称: NDUFB11 Rabbit Monoclonal Antibody
产品货号: AMRe85836



产品概述 (Summary)

产品名称 (Production Name)	NDUFB11 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFB11
别名 (Alternative Names)	ESSS; Np15; P17.3; NP17.3; CI-ESSS; LSDMCA3
基因 ID (Gene ID)	54539.0
蛋白 ID (SwissProt ID)	Q9NX14.Recombinant protein of human NDUFB11

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not


产品名称: NDUFB11 Rabbit Monoclonal Antibody
产品货号: AMRe85836




to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of NDUFB11 in rat heart, rat skeletal muscle lysates using NDUFB11 antibody.




Western blot analysis of NDUFB11 in LnCap, CaCo2 lysates using NDUFB11 antibody



Western blot analysis of NDUFB11 in HeLa, 3T3, PC-3 lysates using NDUFB11 antibody.

产品名称: NDUFB11 Rabbit Monoclonal Antibody
产品货号: AMRe85836



Immunohistochemistry analysis of paraffin-embedded mouse colon using NDUFB11 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .