

产品名称: NEK4 Rabbit Monoclonal Antibody
产品货号: AMRe85845



产品概述 (Summary)

产品名称 (Production Name)	NEK4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEK4
别名 (Alternative Names)	NRK2; STK2; pp12301
基因 ID (Gene ID)	6787.0
蛋白 ID (SwissProt ID)	P51957.A synthetic peptide of human NEK4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 95 kDa; Observed MW: 95 kDa

研究背景 (Background)

Protein kinase that seems to act exclusively upon threonine residues . Required for normal entry into proliferative arrest


产品名称: NEK4 Rabbit Monoclonal Antibody
产品货号: AMRe85845



after a limited number of cell divisions, also called replicative senescence. Required for normal cell cycle arrest in response to double-stranded DNA damage.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NEK4 in Jurkat, HepG2, HeLa lysates using NEK4 antibody.

注意事项 (Note)

For research use only .