

产品名称: Neurogranin Rabbit Monoclonal Antibody
产品货号: AMRe85852



产品概述 (Summary)

产品名称 (Production Name)	Neurogranin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Neurogranin
别名 (Alternative Names)	neurogranin (protein kinase C substrate; RC3); RC3; hng
基因 ID (Gene ID)	4900.0
蛋白 ID (SwissProt ID)	Q92686.Recombinant protein of human Neurogranin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 8 kDa; Observed MW: 14-16 kDa

研究背景 (Background)

Acts as a 'third messenger' substrate of protein kinase C-mediated molecular cascades during synaptic development and


产品名称: Neurogranin Rabbit Monoclonal Antibody
产品货号: AMRe85852



remodeling.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Neurogranin in U87-MG, PC-3 lysates using Neurogranin antibody.

注意事项 (Note)

For research use only .