

产品名称: SLC12A1 Rabbit Monoclonal Antibody
产品货号: AMRe85861



产品概述 (Summary)

产品名称 (Production Name)	SLC12A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC12A1
别名 (Alternative Names)	BSC1; NKCC2
基因 ID (Gene ID)	6557.0
蛋白 ID (SwissProt ID)	Q13621.Recombinant protein of human NKCC2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 121 kDa; Observed MW: 150 kDa

研究背景 (Background)

Electrically silent transporter system. Mediates sodium and chloride reabsorption. Plays a vital role in the regulation of ionic


产品名称: SLC12A1 Rabbit Monoclonal Antibody
产品货号: AMRe85861



balance and cell volume.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SLC12A1/NKCC2 in HeLa lysates using SLC12A1 antibody.

注意事项 (Note)

For research use only .