

产品名称: NONO Rabbit Monoclonal Antibody
产品货号: AMRe85862



产品概述 (Summary)

产品名称 (Production Name)	NONO Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NONO
别名 (Alternative Names)	P54; NMT55; NRB54; P54NRB
基因 ID (Gene ID)	4841.0
蛋白 ID (SwissProt ID)	Q15233.A synthetic peptide of human NMT55/p54nrb

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

研究背景 (Background)

DNA- and RNA binding protein, involved in several nuclear processes. Binds the conventional octamer sequence in double


产品名称: NONO Rabbit Monoclonal Antibody
产品货号: AMRe85862



stranded DNA. Also binds single-stranded DNA and RNA at a site independent of the duplex site (By similarity). Involved in pre-mRNA splicing, probably as a heterodimer with SFPQ.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NMT55/p54nrb in 293T, MOLT4 lysates using NONO antibody.

注意事项 (Note)

For research use only .