

产品名称: THO Complex Subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe85873



产品概述 (Summary)

产品名称 (Production Name)	THO Complex Subunit 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THO Complex Subunit 1
别名 (Alternative Names)	Death domain containing protein p84N5; HPR1; hTRES84; Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1
基因 ID (Gene ID)	9984.0
蛋白 ID (SwissProt ID)	Q96FV9.A synthetic peptide of human Nuclear Matrix Protein p84

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 84 kDa

产品名称: THO Complex Subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe85873




研究背景 (Background)

Regulates transcriptional elongation of a subset of genes. Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of Nuclear Matrix Protein p84 in Jurkat, HeLa, COS7 lysates using THO Complex Subunit 1 antibody.



Immunocytochemistry analysis of Nuclear Matrix Protein p84(green) in HeLa using Nuclear Matrix Protein p84 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .