

产品名称: OCT1 Rabbit Monoclonal Antibody
产品货号: AMRe85883



产品概述 (Summary)

产品名称 (Production Name)	OCT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OCT1 POU2F1; OCT1; OTF1; POU domain; class 2; transcription factor 1; NF-A1;
别名 (Alternative Names)	Octamer-binding protein 1; Oct-1; Octamer-binding transcription factor 1; OTF-1; POU2F2; OCT2; OTF2; POU domain; class 2; transcription factor 2; Lymphoid-restricted im
基因 ID (Gene ID)	5451.0
蛋白 ID (SwissProt ID)	P14859.Recombinant protein of human Oct-1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 90 kDa

产品名称: OCT1 Rabbit Monoclonal Antibody
产品货号: AMRe85883




研究背景 (Background)

Transcription factor that binds to the octamer motif (5'-ATTGTCAT-3') and activates the promoters of the genes for some small nuclear RNAs (snRNA) and of genes such as those for histone H2B and immunoglobulins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Oct1 in HeLa, 293 lysates using Oct1 antibody.

注意事项 (Note)

For research use only .