

产品名称: ORP150 Rabbit Monoclonal Antibody
产品货号: AMRe85887



产品概述 (Summary)

产品名称 (Production Name)	ORP150 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ORP150
别名 (Alternative Names)	Grp170; HSP12A; Hyou1; Orp150
基因 ID (Gene ID)	10525.0
蛋白 ID (SwissProt ID)	Q9Y4L1.A synthetic peptide of human ORP150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 111 kDa; Observed MW: 150 kDa

研究背景 (Background)

Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular


产品名称: ORP150 Rabbit Monoclonal Antibody
产品货号: AMRe85887



chaperone and participate in protein folding.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ORP150 in 293T, MCF-7 lysates using ORP150 antibody.

注意事项 (Note)

For research use only .