

产品名称: p40 Rabbit Monoclonal Antibody
产品货号: AMRe85895



产品概述 (Summary)

产品名称 (Production Name)	p40 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	p40
别名 (Alternative Names)	p40; RAB9p40; RABEPK; RGD1310612; RP11 65N13.1; RP23-446N16.2
基因 ID (Gene ID)	10244.0
蛋白 ID (SwissProt ID)	Q7Z6M1.Recombinant protein of human p40

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)


Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

产品名称: p40 Rabbit Monoclonal Antibody
产品货号: AMRe85895



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of p40 in K562, HeLa lysates using p40 antibody.

注意事项 (Note)

For research use only .