

产品名称: p53 Rabbit Monoclonal Antibody
产品货号: AMRe85898



产品概述 (Summary)

产品名称 (Production Name)	p53 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	p53
别名 (Alternative Names)	bbl; bfy; bhy; p44; p53; Tp53
基因 ID (Gene ID)	22059.0
蛋白 ID (SwissProt ID)	P02340.Recombinant protein of mouse p53

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 53 kDa

研究背景 (Background)

Acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological


产品名称: p53 Rabbit Monoclonal Antibody
产品货号: AMRe85898



circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. One of the activated genes is an inhibitor of cyclin-dependent kinases.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of p53 in 3T3 lysates using p53 antibody.

注意事项 (Note)

For research use only .