

产品名称: AKAP14 Rabbit Monoclonal Antibody
产品货号: AMRe86173



产品概述 (Summary)

产品名称 (Production Name)	AKAP14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AKAP14
别名 (Alternative Names)	AKAP28; A-kinase anchor protein 14; AKAP-14; A-kinase anchor protein 28 kDa; AKAP 28; Protein kinase A-anchoring protein 14; PRKA14;AKAP14
基因 ID (Gene ID)	158798.0
蛋白 ID (SwissProt ID)	Q86UN6.KLH-conjugated synthetic peptide encompassing a sequence within human AKAP14. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	22 kD

产品名称: AKAP14 Rabbit Monoclonal Antibody
产品货号: AMRe86173




研究背景 (Background)

Recombinant rabbit monoclonal antibody to AKAP14

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AKAP14 expression in HEK293T (A), HepG2 (B), Jurkat (C), mouse brain (D), mouse liver (E), rat brain (F), rat liver (G) whole cell lysates. (Predicted band size: 22 kD; Observed band size: 22 kD)

注意事项 (Note)

For research use only .