

产品名称: XPO5 Rabbit Monoclonal Antibody
产品货号: AMRe86174



产品概述 (Summary)

产品名称 (Production Name)	XPO5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	XPO5
别名 (Alternative Names)	KIAA1291; RANBP21; Exportin-5; Exp5; Ran-binding protein 21;Exportin 5
基因 ID (Gene ID)	57510.0
蛋白 ID (SwissProt ID)	Q9HAV4.KLH-conjugated synthetic peptide encompassing a sequence within human Exportin 5. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	136 kD

研究背景 (Background)


产品名称: XPO5 Rabbit Monoclonal Antibody
产品货号: AMRe86174



Recombinant rabbit monoclonal antibody to Exportin 5

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Exportin 5 expression in HEK293T (A), MCF7 (B), Jurkat (C), mouse brain (D), rat brain (E) whole cell lysates. (Predicted band size: 136 kD; Observed band size: 136 kD)

注意事项 (Note)

For research use only .