

产品名称: GTF2F2 Rabbit Monoclonal Antibody
产品货号: AMRe86177



产品概述 (Summary)

产品名称 (Production Name)	GTF2F2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2F2 RAP30; General transcription factor IIF subunit 2; ATP-dependent helicase
别名 (Alternative Names)	GTF2F2; General transcription factor IIF 30 kDa subunit; Transcription initiation factor IIF subunit beta; TFIIF-beta; Transcription initiation factor RAP30; TFIIF RAP 30
基因 ID (Gene ID)	2963.0
蛋白 ID (SwissProt ID)	P13984.KLH-conjugated synthetic peptide encompassing a sequence within human TFIIF RAP 30. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
----------------------	---------------------------------

产品名称: GTF2F2 Rabbit Monoclonal Antibody
产品货号: AMRe86177




蛋白分子量 (Molecular Weight) 28 kD

研究背景 (Background)

Recombinant rabbit monoclonal antibody to TFIIF RAP 30

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TFIIF RAP 30 expression in A549 (A), HeLa (B), K562 (C) whole cell lysates. (Predicted band size: 28 kD; Observed band size: 28 kD)

注意事项 (Note)

For research use only .