

产品名称: KLRG1 Rabbit Monoclonal Antibody
产品货号: AMRe86182



产品概述 (Summary)

产品名称 (Production Name)	KLRG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLRG1 CLEC15A; MAFA; MAFAL; Killer cell lectin-like receptor subfamily G member 1;
别名 (Alternative Names)	C-type lectin domain family 15 member A; ITIM-containing receptor MAFA-L; MAFA-like receptor; Mast cell function-associated antigen; KLRG1
基因 ID (Gene ID)	10219.0
蛋白 ID (SwissProt ID)	Q96E93.KLH-conjugated synthetic peptide encompassing a sequence within human KLRG1. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	38 kD

产品名称: KLRG1 Rabbit Monoclonal Antibody
产品货号: AMRe86182




研究背景 (Background)

Recombinant rabbit monoclonal antibody to KLRG1

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of KLRG1 expression in HEK293T (A), A549 (B), PC3 (C) whole cell lysates. (Predicted band size: 21 kD;
Observed band size: 38 kD)

注意事项 (Note)

For research use only .