

产品名称: ALDH16A1 Rabbit Monoclonal Antibody
产品货号: AMRe86187



产品概述 (Summary)

产品名称 (Production Name)	ALDH16A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALDH16A1
别名 (Alternative Names)	Aldehyde dehydrogenase family 16 member A1;ALDH16A1
基因 ID (Gene ID)	126133.0
蛋白 ID (SwissProt ID)	Q8IZ83.KLH-conjugated synthetic peptide encompassing a sequence within human ALDH16A1. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	80 kD

研究背景 (Background)


产品名称: ALDH16A1 Rabbit Monoclonal Antibody
产品货号: AMRe86187



Recombinant rabbit monoclonal antibody to ALDH16A1

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ALDH16A1 expression in HEK293T (A), HepG2 (B), THP1 (C), mouse liver (D), mouse kidney (E), rat liver (F), rat kidney (G) whole cell lysates. (Predicted band size: 85 kD; Observed band size: 80 kD)

注意事项 (Note)

For research use only .