

产品名称: AGO3 Rabbit Monoclonal Antibody
产品货号: AMRe86192



产品概述 (Summary)

产品名称 (Production Name)	AGO3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGO3 EIF2C3; Protein argonaute-3; Argonaute3; hAgo3; Argonaute RISC catalytic component 3; Eukaryotic translation initiation factor 2C 3; eIF-2C 3; eIF2C3; eIF2C3
别名 (Alternative Names)	
基因 ID (Gene ID)	192669.0
蛋白 ID (SwissProt ID)	Q9H9G7.KLH-conjugated synthetic peptide encompassing a sequence within human eIF2C3. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	97 kD

产品名称: AGO3 Rabbit Monoclonal Antibody
产品货号: AMRe86192




研究背景 (Background)

Recombinant rabbit monoclonal antibody to eIF2C3

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of eIF2C3 expression in HEK293T (A), mouse liver (B) whole cell lysates. (Predicted band size: 97; 71 kD;
Observed band size: 97 kD)

注意事项 (Note)

For research use only .