

产品名称: MYH7B Rabbit Monoclonal Antibody
产品货号: AMRe86197



产品概述 (Summary)

产品名称 (Production Name)	MYH7B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH7B
别名 (Alternative Names)	KIAA1512; Myosin-7B; Antigen MLAA-21; Myosin cardiac muscle beta chain; Myosin heavy chain 7B cardiac muscle beta isoform; Slow A MYH14; MYH7B
基因 ID (Gene ID)	57644.0
蛋白 ID (SwissProt ID)	A7E2Y1.KLH-conjugated synthetic peptide encompassing a sequence within human MYH7B. The exact sequence is proprietary.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	221 kD

产品名称: MYH7B Rabbit Monoclonal Antibody
产品货号: AMRe86197




研究背景 (Background)

Recombinant rabbit monoclonal antibody to MYH7B

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MYH7B expression in mouse liver (A) whole cell lysates. (Predicted band size: 225; 20 kD; Observed band size: 221 kD)

注意事项 (Note)

For research use only .