

产品名称: Pentraxin 3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe86202



产品概述 (Summary)

产品名称 (Production Name)	Pentraxin 3/PTX3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Pentraxin 3/PTX3
别名 (Alternative Names)	TSG-14; TNFAIP5
基因 ID (Gene ID)	5806
蛋白 ID (SwissProt ID)	P26022.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:42 kDa; Observed MW:42 kDa

研究背景 (Background)


产品名称: Pentraxin 3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe86202



This gene encodes a member of the pentraxin protein family. The expression of this protein is induced by inflammatory cytokines in response to inflammatory stimuli in several mesenchymal and epithelial cell types, particularly endothelial cells and mononuclear phagocytes. The protein promotes fibrocyte differentiation and is involved in regulating inflammation and complement activation. It also plays a role in angiogenesis and tissue remodeling. The protein serves as a biomarker for several inflammatory conditions. [provided by RefSeq, Jun 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from 3T3 cells using AMRe86202 at 1:1000.

注意事项 (Note)

For research use only .