

产品名称: LINE 1 ORF1p Rabbit Monoclonal Antibody
产品货号: AMRe86204



产品概述 (Summary)

产品名称 (Production Name)	LINE 1 ORF1p Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LINE 1 ORF1p
别名 (Alternative Names)	L1-ORF1p; LINE-1 retrotransposable element ORF1 protein
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P11260.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:41 kDa; Observed MW:41 kDa

产品名称: LINE 1 ORF1p Rabbit Monoclonal Antibody
产品货号: AMRe86204




研究背景 (Background)

Nucleic acid-binding protein which is essential for retrotransposition of LINE-1 elements in the genome. Functions as a nucleic acid chaperone binding its own transcript and therefore preferentially mobilizing the transcript from which they are encoded.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from F9 cells using AMRe86204 at 1:1000.

注意事项 (Note)

For research use only .