

产品名称: Retinoic Acid Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe86220



产品概述 (Summary)

产品名称 (Production Name)	Retinoic Acid Receptor alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Retinoic Acid Receptor alpha
别名 (Alternative Names)	RAR; NR1B1; RARalpha
基因 ID (Gene ID)	5914
蛋白 ID (SwissProt ID)	P10276.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:51 kDa; Observed MW:55 kDa


产品名称: Retinoic Acid Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe86220

研究背景 (Background)

This gene represents a nuclear retinoic acid receptor. The encoded protein, retinoic acid receptor alpha, regulates transcription in a ligand-dependent manner. This gene has been implicated in regulation of development, differentiation, apoptosis, granulopoiesis, and transcription of clock genes. Translocations between this locus and several other loci have been associated with acute promyelocytic leukemia. Alternatively spliced transcript variants have been found for this locus. [provided by RefSeq, Sep 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, 3T3 cells using AMRe86220 at 1:1000.

注意事项 (Note)

For research use only .