

产品名称: C Reactive Protein Rabbit Monoclonal Antibody
产品货号: AMRe86224



产品概述 (Summary)

产品名称 (Production Name)	C Reactive Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C Reactive Protein
别名 (Alternative Names)	PTX1; CRP
基因 ID (Gene ID)	1401
蛋白 ID (SwissProt ID)	P02741.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:25 kDa; Observed MW:25 kDa

研究背景 (Background)


产品名称: C Reactive Protein Rabbit Monoclonal Antibody
产品货号: AMRe86224



The protein encoded by this gene belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during acute phase response to tissue injury, infection, or other inflammatory stimuli. [provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe86224 at 1:1000.

注意事项 (Note)

For research use only .