

产品名称: Kappa Light chain Rabbit Monoclonal Antibody
产品货号: AMRe86226



产品概述 (Summary)

产品名称 (Production Name)	Kappa Light chain Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Kappa Light chain
别名 (Alternative Names)	HCAK1; Ig kappa chain C region; IGKC; IMMUNOGLOBULIN InV; Km
基因 ID (Gene ID)	3514
蛋白 ID (SwissProt ID)	P01834.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:25 kDa; Observed MW:25 kDa

研究背景 (Background)


产品名称: Kappa Light chain Rabbit Monoclonal Antibody
产品货号: AMRe86226



Constant region of immunoglobulin light chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes. In the recognition phase of humoral immunity, the membrane-bound immunoglobulins serve as receptors which, upon binding of a specific antigen, trigger the clonal expansion and differentiation of B lymphocytes into immunoglobulins-secreting plasma cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe86226 at 1:1000.

注意事项 (Note)

For research use only .