

产品名称: Formyl-Histone H3 (Lys122) Rabbit Monoclonal Antibody
产品货号: AMRe86229



产品概述 (Summary)

产品名称 (Production Name)	Formyl-Histone H3 (Lys122) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Formyl-Histone H3 (Lys122)
别名 (Alternative Names)	H3/A; H3FA
基因 ID (Gene ID)	8350
蛋白 ID (SwissProt ID)	P68431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:15 kDa; Observed MW:17 kDa


产品名称: Formyl-Histone H3 (Lys122) Rabbit Monoclonal Antibody
产品货号: AMRe86229

研究背景 (Background)

Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, 3T3 cells using AMRe86229 at 1:1000.

注意事项 (Note)

For research use only .