

产品名称: SPINT1 Rabbit Monoclonal Antibody
产品货号: AMRe86237



产品概述 (Summary)

产品名称 (Production Name)	SPINT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPINT1
别名 (Alternative Names)	HAI; HAI1; MANSC2
基因 ID (Gene ID)	6692
蛋白 ID (SwissProt ID)	O43278.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:2000,ICC/IF 1:20-1:50,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:58 kDa; Observed MW:45-60 kDa

产品名称: SPINT1 Rabbit Monoclonal Antibody
产品货号: AMRe86237




研究背景 (Background)

The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse large intestine tissue using AMRe86237 at 1:1000.

注意事项 (Note)

For research use only .