

产品名称: Cytokeratin 18 Rabbit Monoclonal Antibody
产品货号: AMRe86238



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 18 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytokeratin 18
别名 (Alternative Names)	K18; CK18; Krt1-18
基因 ID (Gene ID)	16668
蛋白 ID (SwissProt ID)	P05784.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa; Observed MW:48 kDa

研究背景 (Background)


产品名称: Cytokeratin 18 Rabbit Monoclonal Antibody
产品货号: AMRe86238



When phosphorylated, plays a role in filament reorganization. Involved in the delivery of mutated CFTR to the plasma membrane. Involved in the uptake of thrombin-antithrombin complexes by hepatic cells (By similarity). Together with KRT8, is involved in interleukin-6 (IL-6)-mediated barrier protection.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86238 at 1:20000.

注意事项 (Note)

For research use only .