

产品名称: Hsp105 Rabbit Monoclonal Antibody
产品货号: AMRe86261



产品概述 (Summary)

产品名称 (Production Name)	Hsp105 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp105
别名 (Alternative Names)	HSP105; HSP105A; HSP105B; NY-CO-25
基因 ID (Gene ID)	10808
蛋白 ID (SwissProt ID)	Q92598.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:20-1:50, FC 1:20-1:50, IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:97 kDa; Observed MW:105 kDa

研究背景 (Background)


产品名称: Hsp105 Rabbit Monoclonal Antibody
产品货号: AMRe86261




This gene encodes a member of the heat shock protein 70 family of proteins. The encoded protein functions as a nucleotide exchange factor for the molecular chaperone heat shock cognate 71 kDa protein (Hsc70). In addition, this protein plays a distinct but related role as a holdase that inhibits the aggregation of misfolded proteins, including the cystic fibrosis transmembrane conductance regulator (CFTR) protein. Elevated expression of this protein has been observed in numerous human cancers. [provided by RefSeq, Mar 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86261 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling Hsp105 with AMRe86261.

注意事项 (Note)

For research use only .