

产品名称: CTNNBL1 Rabbit Monoclonal Antibody
产品货号: AMRe86266



产品概述 (Summary)

产品名称 (Production Name)	CTNNBL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CTNNBL1
别名 (Alternative Names)	NAP; P14L; PP8304; C20orf33; dJ633O20.1
基因 ID (Gene ID)	56259
蛋白 ID (SwissProt ID)	Q8WYA6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:20-1:50,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:65 kDa; Observed MW:65 kDa

研究背景 (Background)


产品名称: CTNNBL1 Rabbit Monoclonal Antibody
产品货号: AMRe86266




The protein encoded by this gene is a component of the pre-mRNA-processing factor 19-cell division cycle 5-like (PRP19-CDC5L) protein complex, which activates pre-mRNA splicing and is an integral part of the spliceosome. The encoded protein is also a nuclear localization sequence binding protein, and binds to activation-induced deaminase and is important for antibody diversification. This gene may also be associated with the development of obesity. Alternative splicing results in multiple transcript variants. A pseudogene of this gene has been defined on the X chromosome. [provided by RefSeq, Jul 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86266 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling CTNNBL1 with AMRe86266.

注意事项 (Note)

For research use only .